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
FAULT - read any information printed by this symbol.

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2.0 Instalación

2.1 Unpacking and Handling

When unpacking or moving the product always lift it by its base. Do not use the door or any other protruding parts. The product contains a transformer and is heavy: use two or more people to carry the product.

Remove any packing material from the door great and from inside the product chamber. Close the door with care to prevent any damage to insulation surfaces.

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
If the product is to be used in a process which could liberate hazardous gases, then it should be installed together with a suitable fume extraction system.

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2.3 Fitting the Heating Elements

Please refer to the Element Installation and Replacement in section 6.0. Wear eye protection when handling the heating elements.

The molybdenum disilicide elements are **EXCEPTIONALLY FRAGILE** and are packed separately, together with other items shown in the list.

Separately Packed Items	1700 °C Models
Elements	4
Element Clamps	8
Element Clips	8
Braids	1 set of 5
Separators/ Block	4
Chimney Unit	1

The installation of these elements is described in section 6.9

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These models are designed only for single phase electrical supplies, or for two live phases and neutral of a 3-phase supply. The product must be connected only to the type and voltage of supply for which it was ordered.

Check that the supply voltage is compatible with the voltage on the label, and that the current capacity is sufficient for the current on the label, before connection to the electrical supply. A table of the most common ratings is given in section 9.0 towards the back of this manual.

Normally a supply cable is not fitted, and connection is to be made to terminal blocks accessed by removal of the left-hand end side cover (see section 6.4). The electrical supply cable must be properly connected and fitted with an appropriate strain relief at the entry to the case.

The electrical supply cable should be wired to an isolator or fitted with a line plug. The isolator must be within easy reach of the operator, or the line plug must be easily removable and on a lead not exceeding 3 m length. The line plug must allow the protective earth (ground) to be connected before the supply conductors and to be disconnected after them. The supply point must be marked as the disconnecting device for the furnace.

The electrical supply **MUST** incorporate an earth (ground).

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2-phase	L1	Black	to phase 1
	L2	Black	to phase 2
	N	Light Blue	to neutral
	PE	Green/ Yellow	to earth (ground)

2.5 Power Supply Notes

Each model is manufactured for one of the following voltage ranges:

- 208 V
- 220 to 240 V
- 380 to 415 V (not HTF 17/5 or HTF 18/4)

It is not possible to modify and use a furnace manufactured for the 380-415 range on either of the other two voltage ranges listed above: there are too many component differences.

It is possible to modify a furnace manufactured for either 208 and 220-240, provided the thyristor unit is replaced. A transformer tapping must also be altered. See sections 6.5 and 6.11.

It is possible to alter the voltage within any of the ranges above, by reconnecting the incoming cable to the transformer to the appropriate primary tapping. The tapings are 208 - 220 - 230 - 240 V, or 380 - 400 - 415 V, and are labelled on the primary side of the transformer - see section 6.11. It is also necessary to adjust the setting of the thyristor stack: see section 6.5.

Examples:

- to alter a furnace made for a 240 V supply to 208 V: replace thyristor stack, move a cable to the 208 V transformer tapping, and adjust the thyristor stack.
- to change from a supply of 220 V to 230 V: move a cable to the 230 V, and adjust the thyristor stack.

Please contract Carbolite Gero Service for guidance and assistance if the power supply shown on the rating label does not match the power supply available.

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4.0 Funcionamiento

4.1 Operating Cycle

The product is fitted with an instrument switch. The switch cuts off power to the controllers and contactor.

Connect the product to the electrical supply. The cooling fans will operate.

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The heat light glows brightly at first, more dimly as the product temperature approaches a program setpoint.

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To switch the product off, turn the instrument switch off. The case cooling fans continue to operate. Leave the fans on until the product cools to below 300 °C. If the product is to be left unattended, isolate it from the electrical supply.

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4.3 Loading The Furnace

When heating large objects, in particular poor conductors, avoid shielding the thermocouple from the heating elements. Also ensure that nothing is within 15 mm of the elements.

The thermocouple is intended to sense the temperature near the heating element. However if a large object is placed in the chamber it may record the average temperature of the object and the elements, which can lead to over heating of the elements. Allow large objects to gain heat at a lower temperature and then reset the controller to a temperature close to the desired maximum.

Take care that nothing hits the elements when loading and unloading.

4.4 Opening the Door

Take great care when loading or unloading the furnace chamber. See section 4.2. Before removing a hot object from the furnace, ensure a suitable surface is available on which to put it.

Do not open the door at high temperatures. If possible, do not open it above 200 °C. If it is necessary to load or unload work at elevated temperatures, keep the door open for as short a period as possible. The insulation cools quickly and may crack through thermal shock.

The heating elements and the insulation are very susceptible to mechanical shock. At all times operate the door gently to avoid mechanical shock.

Opening the door when the furnace is very hot can cause hot air to be sucked by the case cooling fans towards the thermal cut out and activate it (see section 4.11). Avoid excessive door opening and shut it as soon as possible.

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4.6 Atmospheres & Corrosive Materials

In oxidising atmospheres, metal oxides react with the silica layer on the surface of the elements and may lead to premature failure. Protect the elements from splashes of molten metal and dust when loading the furnace and from fumes developed when melting, especially from fluxes. Also avoid compounds with a high alkali content.

The furnace is designed for use up to its maximum temperature in oxidising atmospheres, but can be operated successfully in neutral or carburising atmospheres.

It may be used with nitrogen, argon or helium atmospheres to 1600 °C. Reducing atmospheres are not recommended. Sulphur dioxide is not harmful in low concentrations. However chlorine and fluorine attacks oxidised elements and should be avoided. A harmful gas is produced and collects in poorly ventilated spaces.

The furnace is not recommended for burning off carbonaceous materials. Other Carbolite Gero furnaces are available for this application.

When an optional gas inlet is fitted there is a label near the inlet saying "INERT GAS ONLY".

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Periodically operate the furnace in air at 1500 °C for an hour to ensure the protective glaze on the elements surface is maintained.

The hearth can be protected from abrasion, if required, by the provision of a secondary plate. This consumable item can either be in lightweight ceramic fibre or in dense, hard-wearing alumina refractory.

4.7 Pesting

The furnace is fitted with molybdenum disilicide elements and is therefore susceptible to peesting. Pesting can be seen as a yellowish crystal growth on the surface of the element. It occurs when the element surface temperature is around 450 °C and is caused by the air oxidising with the molybdenum in the element compound. Under normal operating conditions (above 800 °C) peesting is not a problem and can often be seen occurring naturally on the element terminals.

operating the furnace for prolonged periods at temperature below 600 °C will cause peesting to occur on the element surface and will eventually damage the elements.

Heating on a slow ramp rate or dwelling at temperatures below 600 °C for short periods will not cause a problem providing the furnace is operated above 1400 °C at regular intervals to recondition the element surfaces.

Large amounts of peesting occurring on the element terminals must be removed in accordance with good health and safety practices.

4.8 Explosive Materials

The product must not be used to heat materials which could explode, or which could emit gases that could form explosive mixtures.

4.9 Notes on Temperature Control

This product is designed for heating and cooling at controlled rates. The heating and cooling rates can be set in the temperature programmer. The programmer enables the furnace to heat or cool at slower rates as desired and variable "hold" (dwell) periods can be programmed as required.

The programmer is used in conjunction with a phase angle thyristor power controller, which incorporates a current limit potentiometer pre-set by Carbolite Gero, but which will require adjustment in the event of change of supply voltage.

The elements are connected in series across the low voltage output of a transformer housed in the furnace case. Molybdenum disilicide elements do not age, if an element fails it is not necessary to replace the complete set.

4.10 Thermocouple Warnings

(1) The output from 1700-1800 °C thermocouples when used regularly at temperatures greater than 1650 °C can deteriorate and decrease with age faster than if used at temperatures below 1650 °C; this will cause the furnace to operate at temperatures higher than indicated.

Operators are advised to periodically check the thermocouple output, either by a calibration test or by comparing the output with a new reference thermocouple which has been subjected to high temperatures for a minimum length of time.

Failure to check the thermocouple regularly may result in overheating of the work and the furnace, with consequential damage to both.

(2) The thermocouples fitted to these models give very low outputs below approximately 600 °C and do not give accurate readings at low temperatures. They may show a negative temperature when the furnace is started from cold.

This furnace is not intended to operate with a setpoint below 600 °C.

4.11 Thermal Cutouts

There is a thermal cutout on the interior case surface at the centre top near the chimney. In the event of fan failure or any other reason for case over heating, power to the heating elements is cut.

To reset the thermal cutout, first isolate the furnace from the electrical supply. Remove the top panel and press upwards the button on the cutout device. TAKE CARE - the chimney may be hot.

There is a further thermal cutout in the element circuit transformer which cuts power if the transformer core overheats. This is self-resetting.

If either thermal cutout is activated, then a fault light on the control panel is illuminated.

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The product's elements are very susceptible to mechanical shock. Take great care when loading or unloading the chamber. If it is necessary to load or unload work at elevated temperatures, keep the door open for as short a period as possible. The insulation cools quickly and may crack though thermal shock.

For improved insulation and element life it is recommended to heat and cool at a slow ramp rate e.g. 5 °C/minute, and to avoid opening the door at high temperatures.

On first installing the elements and on subsequent element replacement, operate the product at 1500 °C for an hour to create a protective glaze on the element surface.

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5.3.1 Other Electrical Components

Regular visual or electrical checks should be made on the condition of the electrical supply cable. All internal fuses and visible internal cables should be inspected periodically.

5.3.2 Element Glaze

Depending on the use, heating elements may lose their glaze and gain a rough appearance. They should be checked from time to time. If the glaze has disappeared, it may be restored by heating the furnace up to 1500 °C, without load, for approximately 2 hours.

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5.6 Recommended Spare Parts and Spare Parts Kit

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Each kit consists of one thermocouple, one sheath, one power thyristor, one door insulation piece and a set of elements, clips and braids, element clamps and insulators. Individual spare parts are also available.

When ordering spare parts please quote the model details as requested above.

5.7 Power Adjustment (Controller)

The product's controller incorporates a power limit parameter OP.Hi, which is usually inaccessible to the operator.

The correct setting for the power limit depends on the supply voltage, table given in section 9.0

5.8 Power Adjustment (Thyristor)

The current-limiting thyristor stacks which control power to the elements are fitted with an adjustable resistor which is factory set to limit the maximum current supplied. In the event of a change of supply voltage, or the fitting of a new thyristor, further adjustment may be required.

The maximum element currents for this model are listed in section 9.0. Please contact Carbolite Gero for further information.

5.9 Low Voltage Compensation

If the supply voltage proves to be routinely below the nominal figure for which the furnace has been set up, there is a "+2%" position on the transformer primary side which may be used instead of the standard position, see section 6.11.

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6.3 Safety Warning - Molybdenum Disilicide Elements

Molybdenum disilicide elements form a glazed surface when heated. Internal stresses can form through heating and cooling which render the glaze fragile. The glaze can sometimes splinter into a shower of sharp particles when handled. Always wear eye protection when handling the elements.

Handle the heating elements with extreme care as they are very fragile. Also, avoid touching the heating surface (the thin part of the element), as the material is susceptible to corrosive damage from skin contact.

6.4 Side Panel Removal

Except where explicitly stated, always disconnect the electrical supply before removing the side panel.

Remove the panel by loosening the four fixing screws (behind plastic caps) at the left-hand end of the furnace; do not remove the screws. Lift the panel about 15 mm and then pull off to the side.

6.5 Thyristor Replacement and Adjustment

Replacement

To replace the thyristor unit, isolate the furnace from the electrical supply and remove the left-hand side cover. See section 6.4. Make a note of all wiring to the thyristor, then disconnect it. Replace the unit and connect the wiring again.


Note that the thyristor unit requires replacement if changing from, or to, a 208 V supply.

If the thyristor unit has been replaced because of a change of voltage, the correct transformer tapping should also be set. See section 6.11 for panel removal instructions.

Adjustment

After any thyristor unit replacement, or any change of voltage or transformer tapping, the potentiometer on the thyristor must be adjusted to give the correct element current. This should be done by a qualified person, as there are dangerous voltages inside the control compartment. It also requires a calibrated non-intrusive clip-on ammeter.

Before connecting the electrical supply, turn the potentiometer on the thyristor fully to the left (anti-clockwise). This sets the output current of the thyristor to 'off'.

Connect the electrical supply with the side cover off. TAKE CARE! Set the furnace temperature to maximum. Allow the furnace to start heating up.

Measure the current through the element circuit. This is carried out with the clip-on meter around one pair of thick cables on the left-hand side of the transformer (as seen when facing the control side of the furnace).

Adjust the potentiometer on the face of the thyristor unit. Adjust it slowly to the right (clockwise) to increase the current, pausing to allow the time for response at the meter. Keep adjusting to obtain an ammeter reading of between 149 to 150 A for the HTF 1700 or 139 to 140 A for the HTF 1800. This adjustment should be set up within the first 5 minutes of heat up from ambient and should finally be checked when the furnace temperature is approximately 100 °C below its maximum temperature. Make further adjustments if necessary at this temperature.

Disconnect the electrical supply to ensure safe replacement of the side panel. See section 6.1

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6.7 Fuse Replacement

Access to internal fuses is by removal of the furnace side cover (see section 6.4). See section 9.0 for details of fuses fitted.

6.8 Thermocouple Replacement

Disconnect the product from the supply and remove the product's element access panel.

Make a note of the thermocouple connections. The negative leg of the thermocouple is marked blue. The "compensating" cable for 1700 & 1800 °C thermocouples is plain copper.

Disconnect the thermocouple from its terminal block.

Undo the screw to release the thermocouple sheath; withdraw the sheath and shake out any fragments of thermocouple.

Re-assemble with a new thermocouple observing the colour coding. Ensure that the thermocouple is not twisted as it is being inserted and that the metal tag is bent back, or the screw inserted, to grip the sheath.

6.9 Element Installation and Replacement

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See section 6.3 - Safety Warning - molybdenum disilicide.

Molybdenum disilicide elements form a glazed surface when heated. Internal stresses can form through heating and cooling which render the glaze fragile. The glaze can sometimes splinter into a shower of sharp particles when handled. Always wear eye protection when handling the elements.

Handle the heating elements with extreme care as they are very fragile. Also, avoid touching the heating surface (the thin part of the element), as the material is susceptible to corrosive damage from skin contact.

Initial Installation:

The elements, clips and braids are packed separately. Handle them carefully while unpacking.

Fit the insulation blocks to the elements. Locate the clamps over the elements and tighten carefully. The length of element above the fixing clamps should be as follows:

1700 °C models: 42 mm

1800 °C models: 45 mm

Lower the elements into position and connect the braids according to the scheme shown below using the clip tool provided. The braids must be held tightly to the element as the clips are fitted:

good contact is essential; poor contact can lead to sparking and destruction of the top of the element.

Ensure that the elements are correctly placed: the thin part of the element, and the tapered section, should ideally be within the heating chamber. The element should not touch the bottom of the chamber.

Important -There should be a minimum clearance of 10 mm below the bottom of the element, and a minimum clearance of 15 mm from the sides of the chamber. The elements should be parallel with the side walls.

Replacements:

Read the section above on initial installation.

Remove the aluminium braids and clips using the clip tool provided with the furnace. Lift out the old elements and the insulation blocks: handle the insulation pieces with care as they are fragile.


Carefully loosen and remove the clamps.

Prepare and install the new elements as explained above.

Heat up the furnace to a moderate temperature and ensure that the furnace is controlling properly, in case the previous element failure resulted from a fault in the control system.

After installing new elements, run the furnace at 1500 °C for an hour. This creates a protective glaze on the element surface.


Element Fitting Layout

Use of Clip Tool


The tool comprises two levers

1. Clip tool
2. Apply hand pressure

Key	
A	Braid
B	Roof Insulation
C	Element block
D	Element
E	Element Clamp
*	See text

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After any replacement of insulation material, run the furnace at 1500 °C to burn off volatile matter. Do this in a well ventilated area. Try to ensure that there is some chamber ventilation, but not too much as this could result in cracked insulation.

6.11 Transformer Tappings


A change of transformer tapping is required if the electrical supply voltage differs from the present setting.

To access the transformer tappings, first disconnect the furnace from the electrical supply and then remove the left-hand side cover (see section 6.4). Positions below are given from a viewing position at the left end of the furnace.

On the left side are the secondary tappings, which should not require change, but can be checked against the following data:

HTF 17/5	27 V
HTF 17/10	38 V
HTF 18/4	31.1 V
HTF 18/8	43.8 V

On the right side is a terminal strip with the following connections (208 to 240 V version shown).

In the case of the 380-415 V version, there are three tappings marked 380, 400 and 415.

The thermal trip connection forms part of the circuit to the contactor coil.

The L and N connection may alternatively be L1 and L2 for a live-to-live supply. The L side should be connected to the correct tapping to match the electrical supply voltage.

The link wire should be in position as shown unless the electrical supply voltage is always low, in which case there is the option of moving one end to the +2% terminal.

Important - Changing the transformer primary tapping (including the +2%) requires adjustment of the thyristor current - see section 6.5.


The terminal between the link connection is not used.

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
7.0 Fault Analysis

A. Furnace Does Not Heat Up					
1.	The HEAT light(s) are ON.	▶	An ohm meter applied to the element circuit shows an open circuit.	▶	A heating element has failed.
2.	The HEAT light(s) are OFF.	▶	The controller shows a very high temperature or a code such as S.br.	▶	The thermocouple has broken or has a wiring fault.
		▶	The controller shows a low temperature.	▶	The door switch may be faulty or need adjustment.
				▶	The thyristor fuse may have blown.
				▶	The power thyristor could be failing to switch on due to internal failure, faulty wiring from the controller, or fault controller.
		▶	There are no lights glowing on the controller.	▶	Check the supply fuses and any fuses in the product control compartment.
				▶	The controller may be faulty or not receiving a supply due to a faulty switch or a wiring fault.
		▶	The fault light on the control panel is on	▶	The thermal cutout on the case or in the transformer has activated.

B. Product Overheats					
1.	The HEAT light goes OFF with the instrument switch.	▶	The controller shows a very high temperature.	▶	The controller is faulty.
		▶	The controller shows a low temperature.	▶	The thermocouple may have been shorted out or may have been moved out of the furnace.
				▶	The thermocouple may be connected the wrong way around.
				▶	The controller may be faulty.
2.	The HEAT light does not go off with the instrument switch and the fault persists when a 2 A control fuse is removed from its fuse-holder.	▶	The power thyristor has failed "ON".	▶	Check for an accidental wiring fault which could have overloaded the thyristor. <i>Isolate the furnace if this fault persists.</i>

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
	Thyristor Fuse		Ferraz Protistor of the rating shown
--	----------------	--	--------------------------------------

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Model	Phases	Volts	Supply Fuse Rating	Thyristor Fuse Rating	Current Limit (element circuit)†
HTF 17/10	1-phase	208	32 A	50 A	150 A
HTF 17/10	1-phase	220-240	32 A	50 A	150 A
HTF 17/10	2-phase + N*	380-415	20 A	30 A	150 A

* 2-phase designs use the line-to-line voltage; neutral is used for the control circuit

† see section 5.8 - the current limit setting at lower temperatures is 146 A.

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* Maximum power as measured, included ancillary components and transformer losses.

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altitude:	not exceeding 2000 m
electrical supply:	fluctuation not exceeding 10%
overvoltage:	category II IEC60364-4-443
pollution:	degree 2

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Copyright © 2024 Carbolite Gero Limited
Carbolite Gero Ltd,